DE ID 213 Instrumenten-Desinfektion

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EN ID 213 Instrument disinfection

ID 213 is an aldehyde-free concentrate for the disinfection and cleaning of general dental instruments (mirrors, probes, tweezers, forceps etc.). An additional field of application is the disinfection and cleaning of alkali- and alcohol-sensitive rotary instruments.

Composition: 100 g ID 213 contain 15 g 3-aminopropyl-dodecyl-1,3-propanediamine, 12,5 g alkyl-benzyl-dimethyl-ammoniumchloride, non-ionic surfactants, complexing agents and accessory agents.

Effective spectrum: Bactericidal, tuberculocidal, fungicidal¹⁾, virucidal activity²⁾ (encapsulated viruses such as vaccinia viruses, incl. HBV, HCV, HIV²⁾ as well as non-encapsulated viruses, such as adeno viruses³⁾, polyoma viruses SV 40³⁾, polio viruses³⁾. **List of VAH/DGHM.** Tested in accordance with EN 13727, EN 13624, EN 14348, EN 14476, EN 14561, EN 14562, EN 14563.

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☞ see bottle imprint.

Danger: May cause damage to organs through prolonged or repeated exposure. Causes severe skin burns and eye damage. Harmful if swallowed. Very toxic to aquatic life. Harmful to aquatic life with long lasting effects. Avoid release to the environment. Wear protective gloves and eye/face protection. If skin irritation or rash occurs: Get medical advice/attention. IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. Store in a well-ventilated place. Keep container tightly closed. Dispose of contents/container to hazardous or special waste collection point.

DÜRR DENTAL (Products) UK Ltd.
14 Linnway Way - Telford Way Industrial Estate
Kettering, Northants NN16 8PS, United Kingdom ☎ +44 1536526740

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DÜRR DENTAL AG